



Product Specification
Lid for sample tube 73.641 & 73.646, length: 8 mm, LD-PE

Page 1



Product description

Order number	65.649
Product description	Lid for sample tube 73.641 & 73.646, natural, length: 8 mm, material: LD-PE, 1,000 piece(s)/bag

Product characteristics

Cap	Skirted push cap
Suitable for	sample tube 73.641 & 73.646

Size

Length of product	8 mm
-------------------	------

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	natural

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Lid for sample tube 73.641 & 73.646, length: 8 mm, LD-PE

Page 2

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Minimum order qty.	10000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	10000
Piece(s) / pallet	420000
Depth of box	350 mm
Width of box	300 mm
Height of box	45 mm
Depth of case	375 mm
Width of case	378 mm
Height of case	250 mm
Case volume	0.0354 cbm
Weight of product	0.0005 kg
Weight of case	5.53 kg
EAN of inner box	4038917083119
EAN case	4038917002073

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-13

This document was prepared by EDP support and is valid without signature.



SARSTEDT