



### **Product description**

Order number	65.718
Product description	Lid for sample tube 73.666, natural, length: 8 mm, material: LD-PE, 1,000 piece(s)/bag

### **Product characteristics**

Cap	Skirted push cap
Closure type	push cap
Suitable for	sample tube 73.666

### **Size**

Length of product	8 mm
-------------------	------

### **Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	natural

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

### **Packaging**

Minimum order qty.	10000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	10000
Piece(s) / pallet	360000
Depth of box	320 mm
Width of box	320 mm
Height of box	50 mm
Depth of case	398 mm
Width of case	398 mm
Height of case	310 mm
Case volume	0.0491 cbm
Weight of product	0.0006 kg
Weight of case	6.97 kg
EAN of inner box	4038917083669
EAN case	4038917002561

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

