



### **Product description**

Order number	65.718
Product description	Lid for sample tube 73.641, natural, length: 8 mm, material: LD-PE

### **Product characteristics**

Cap	Skirted push cap
Type of cap	push cap
Suitable for	sample tube 73.641

### **Size**

Length of product	8 mm
-------------------	------

### **Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	natural

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Product in accordance with MDD or IVDD	IVD product
CE certified	CE
Lot or serial number	yes

### **Packaging**

Minimum order qty.	10000
Piece(s) / inner box	1000
Piece(s) / case	10000
Piece(s) / pallet	360000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0,0475 cbm
Weight of product	0,0006 kg
Weight of case	6,97 kg
EAN of inner box	4038917083669
EAN case	4038917002561

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

