



Product Specification
Insert, suitable for tubes and S-Monovette® Ø 15 mm, transparent

Page 1



Product description

Order number	73.1052
Product description	Insert, suitable for tubes and S-Monovette® Ø 15 mm, material: PS, length: 23.5 mm, transparent, 500 piece(s)/bag

Product characteristics

Graduation	no
Suitable for	tubes and S-Monovette® Ø 15 mm
Base shape	skirted conical base

Size

Diameter of opening	11 mm
Length of product	23.5 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Insert, suitable for tubes and S-Monovette® Ø 15
mm, transparent

Page 2

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	156000
Depth of case	428 mm
Width of case	323 mm
Height of case	130 mm
Case volume	0.018 cbm
Weight of product	0.0008 kg
Weight of case	2.02 kg
EAN of inner box	4038917085151
EAN case	4038917003834

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT