



Product description

| | |
|---------------------|---|
| Order number | 73.1053 |
| Product description | Insert, suitable for tubes and S-Monovette® Ø 16 mm, material: PS, length: 23.5 mm, transparent, 500 piece(s)/bag |

Product characteristics

| | |
|--------------|--------------------------------|
| Graduation | no |
| Suitable for | tubes and S-Monovette® Ø 16 mm |
| Base shape | skirted conical base |

Size

| | |
|---------------------|---------|
| Diameter of opening | 11 mm |
| Length of product | 23.5 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Insert, suitable for tubes and S-Monovette® Ø 16
mm, transparent

Page 2

Material & colours

| | |
|-------------------|------------------|
| Product material | Polystyrene (PS) |
| Colour of product | transparent |

Purity & certification

| | |
|------------------|---------------|
| Product category | not regulated |
| CE certified | no CE mark |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 2000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 2000 |
| Piece(s) / pallet | 156000 |
| Depth of case | 428 mm |
| Width of case | 323 mm |
| Height of case | 130 mm |
| Case volume | 0.018 cbm |
| Weight of product | 0.0009 kg |
| Weight of case | 2.26 kg |
| EAN of inner box | 4038917085168 |
| EAN case | 4038917003841 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

