



### **Product description**

Order number	73.641
Product description	Sample tube, suitable for auto-analyzer, material: PS, 1.5 ml, (LxØ): 24 x 14 mm, transparent, 1,000 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	push cap
Suitable for	auto-analyzer
Base shape	conical base

### **Size**

Max. volume	1.5 ml
Diameter	14 mm
Diameter of opening	8 mm
Height of product	24 mm
Length of product	24 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Sample tube, suitable for auto-analyzer,**  
**transparent**

**Page 2**

**Material & colours**

Product material Polystyrene (PS)  
Colour of product transparent

**Purity & certification**

Product category In-vitro diagnostic, CE  
CE certified CE  
Batched yes

**Packaging**

Minimum order qty. 5000  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 1000  
Piece (s) / outer case 5000  
Piece(s) / pallet 210000  
Depth of case 388 mm  
Width of case 378 mm  
Height of case 250 mm  
Case volume 0.0367 cbm  
Weight of product 0.001 kg  
Weight of case 5.37 kg  
EAN of inner box 4038917083072  
EAN case 4038917002042

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2024-03-04

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**