



Product description

Order number	73.641
Product description	Sample tube, suitable for auto-analyzer, material: PS, 1.5 ml, (LxØ): 24 x 14 mm, transparent, 1,000 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Suitable for	auto-analyzer
Base shape	conical base

Size

Max. volume	1.5 ml
Diameter	14 mm
Diameter of opening	8 mm
Height of product	24 mm
Length of product	24 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 5000
Type of smallest bag
subpackaging
Piece(s) / inner box 1000
Piece (s) / outer case 5000
Piece(s) / pallet 210000
Depth of case 388 mm
Width of case 378 mm
Height of case 250 mm
Case volume 0.0367 cbm
Weight of product 0.001 kg
Weight of case 5.37 kg
EAN of inner box 4038917083072
EAN case 4038917002042

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

