



Product description

Order number	73.646
Product description	Sample tube, suitable for auto-analyzer, material: PS, 3.5 ml, (LxØ): 38 x 14 mm, transparent, 1000 piece(s)/inner box

Product characteristics

Graduation	no
Type of cap	push cap
Suitable for	auto-analyzer
Base shape	conical base

Size

Max. volume	3,5 ml
Diameter	14 mm
Diameter of opening	8 mm
Height of product	38 mm
Length of product	38 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Sample tube, suitable for auto-analyzer,
transparent

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product in accordance with IVD product
MDD or IVDD
CE certified CE
Lot or serial number yes

Packaging

Minimum order qty. 5000
Type of inner packaging inner box
Piece(s) / inner box 1000
Piece(s) / case 5000
Piece(s) / pallet 120000
Depth of case 398 mm
Width of case 398 mm
Height of case 415 mm
Case volume 0,0657 cbm
Weight of product 0,0018 kg
Weight of case 9,62 kg
EAN of inner box 4038917083096
EAN case 4038917002066

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-04-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT