



### **Product description**

Order number	73.646
Product description	Sample tube, suitable for auto-analyzer, material: PS, 3.5 ml, (LxØ): 38 x 14 mm, transparent, 1000 piece(s)/inner box

### **Product characteristics**

Graduation	no
Type of cap	push cap
Suitable for	auto-analyzer
Base shape	conical base

### **Size**

Max. volume	3,5 ml
Diameter	14 mm
Diameter of opening	8 mm
Height of product	38 mm
Length of product	38 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

### **Purity & certification**

Product in accordance with MDD or IVDD	IVD product
CE certified	CE
Lot or serial number	yes

### **Packaging**

Minimum order qty.	5000
Type of inner packaging	inner box
Piece(s) / inner box	1000
Piece(s) / case	5000
Piece(s) / pallet	120000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0,0657 cbm
Weight of product	0,0018 kg
Weight of case	9,62 kg
EAN of inner box	4038917083096
EAN case	4038917002066

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

