



### **Product description**

Order number	73.650.500
Product description	Sample tube, suitable for Coulter Counter, cap enclosed, material: PS, (LxØ): 60 x 34 mm, 250 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	push cap
Suitable for	Coulter Counter
Base shape	flat base

### **Size**

Diameter	34 mm
Diameter of opening	12 mm
Length of product	60 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polystyrene (PS)
Closure material	Low Density Polyethylene (LD-PE)

### **Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

### **Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	500
Piece(s) / pallet	18000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0077 kg
Weight of case	4.51 kg
EAN of inner box	4038917157841
EAN case	4038917157834

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

