



Product description

Order number	73.663
Product description	Sample tube, suitable for Boehringer ES 300® and ES 600®, material: PS, 2 ml, (LxØ): 40 x 12 mm, transparent, 500 piece(s)/bag

Product characteristics

Graduation	no
Type of cap	push cap
Suitable for	Boehringer ES 300® and ES 600®
Base shape	conical base

Size

Max. volume	2 ml
Diameter	12 mm
Diameter of opening	10 mm
Height of product	40 mm
Length of product	40 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	In-vitro diagnostic, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	5000
Type of inner packaging	bag
Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	180000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0,0475 cbm
Weight of product	0,0012 kg
Weight of case	6,72 kg
EAN of inner box	4038917083195
EAN case	4038917002158

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

