



Product Specification

**Liquor/CSF-Monovette® , 10 ml, cap blue, (LxØ):
102 x 15 mm, 1 piece(s)/blister**

Page 1

Product description

| | |
|---------------------|---|
| Order number | 10.252.027 |
| Product description | Liquor/CSF-Monovette® , nominal volume: 10 ml, without preparation, collection and storage of diagnostic specimens, (LxØ): 102 x 15 mm, transparent, graduated, material: PP, Luer screw cap, cap blue, with paper label, "Liquor/Cerebrospinal Fluid", sterile, 1 piece(s)/blister |

Product characteristics

| | |
|-----------------------|--|
| Type of preparation | without preparation |
| Label/ Print | with paper label |
| Colour of print/label | blue |
| Graduation | yes |
| Cap | Luer screw cap |
| Closure type | screw cap |
| Light protection | no |
| Application | collection and storage of diagnostic specimens |
| Example of use | CSF, collection, transport and storage of CSF/liquor |
| Base shape | round base |
| Volume of work | 10 ml |
| Connection | Luer (m) |

Size

| | |
|----------------------------------|--------|
| Sample volume | 10 ml |
| Diameter | 15 mm |
| Length including cap | 117 mm |
| Length excluding cap | 102 mm |
| Length including cap and plunger | 127 mm |

Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | blue |
| Tube material | Polypropylene (PP) |
| Piston material | High Density Polyethylene (HD-PE) |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Liquor/CSF-Monovette® , 10 ml, cap blue, (LxØ):
102 x 15 mm, 1 piece(s)/blister

Page 2

Piston plunger material Polystyrene (PS)

Purity & certification

| | |
|---------------------------|--------------------------------------|
| Satisfies the requirement | IATA, ADR |
| Product category | In-vitro diagnostic (IVD) |
| CE certified | CE - manufacturer's self-declaration |
| Purity standard | sterile |
| Sterilisation | Electron irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | blister |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 13500 |
| Depth of box | 351 mm |
| Width of box | 151 mm |
| Height of box | 171 mm |
| Depth of case | 778 mm |
| Width of case | 368 mm |
| Height of case | 195 mm |
| Case volume | 0.0558 cbm |
| Weight of product | 0.0077 kg |
| Weight of case | 5.46 kg |
| EAN of inner box | 4038917104487 |
| EAN case | 4038917040334 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-13

This document was prepared by EDP support and is valid without signature.



SARSTEDT