



Product Specification
Micro sample tube 100 µl brilliant, 100 µl, screw cap

Page 1



Product description

Order number	42.117
Product description	Micro sample tube 100 µl brilliant, preparation: 1% brilliant cresyl blue solution, nominal volume: 100 µl, (LxØ) with cap: 46 x 10.8 mm, screw cap, cap: blue, conical base, 50 piece(s)/case

Product characteristics

Type of collection	Venous and capillary
Type of preparation	1% brilliant cresyl blue solution
Preparation concentration	1% brilliant cresyl blue solution
Closure type	screw cap
Example of use	Reticulocyte counting
Base shape	conical base
Volume of work	100 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Micro sample tube 100 µl brilliant, 100 µl, screw cap

Page 2

Size

Sample volume	100 µl
Diameter	10.8 mm
Length of product	46 mm
Length including cap	46 mm
Length excluding cap	44 mm

Material & colours

Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	blue

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	117000
Depth of box	190 mm
Width of box	130 mm
Height of box	15 mm
Depth of case	283 mm
Width of case	173 mm
Height of case	135 mm
Case volume	0.0066 cbm
Weight of product	0.0015 kg
Weight of case	0.92 kg
EAN of inner box	4038917080897
EAN case	4038917365130

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

