



**Product Specification**  
**Cuvette and steel balls, suitable for Amelung KC  
coagulometer (macro), transparent**

**Page 1**



**Product description**

Order number	73.834
Product description	Cuvette and steel balls, suitable for Amelung KC coagulometer (macro), material: PS, length: 23 mm, transparent, 1,000 cuvettes/bag and 1,000 steel balls/flask

**Product characteristics**

Graduation	no
Suitable for	Amelung KC coagulometer (macro)
Base shape	flat base

**Size**

Diameter of opening	10 mm
Length of product	23 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Cuvette and steel balls, suitable for Amelung KC**  
**coagulometer (macro), transparent**

**Page 2**

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

**Purity & certification**

Product category	not regulated
CE certified	no CE mark
Batched	yes

**Packaging**

Minimum order qty.	2000
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	216000
Depth of case	308 mm
Width of case	223 mm
Height of case	195 mm
Case volume	0.0134 cbm
Weight of product	0.0008 kg
Weight of case	1.8 kg
EAN of inner box	4038917084321
EAN case	4038917003186

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-14

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**