



Product description

Order number	73.834.001
Product description	Steel ball for Amelung KC coagulometer (macro), grey

Product characteristics

Suitable for	Amelung KC coagulometer (macro)
--------------	---------------------------------

Material & colours

Colour of product	grey
-------------------	------

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	1000
Piece(s) / inner box	1000
Piece (s) / outer case	50000
Piece(s) / pallet	5400000
Depth of case	308 mm
Width of case	223 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Steel ball for Amelung KC coagulometer (macro)

Page 2

Height of case	195 mm
Case volume	0.0134 cbm
Weight of product	0.0001 kg
Weight of case	7.15 kg
EAN of inner box	4038917121682
EAN case	4038917017930

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-24

This document was prepared by EDP support and is valid without signature.



SARSTEDT