



Product description

Order number	73.834
Product description	Cuvette and steel balls, suitable for Amelung KC coagulometer (macro), material: PS, length: 23 mm, transparent, 1,000 cuvettes/bag and 1,000 steel balls/flask

Product characteristics

Graduation	no
Suitable for	Amelung KC coagulometer (macro)
Base shape	flat base

Size

Diameter of opening	10 mm
Length of product	23 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 2000
Piece(s) / inner box 1000
Piece (s) / outer case 2000
Piece(s) / pallet 216000
Depth of case 308 mm
Width of case 223 mm
Height of case 195 mm
Case volume 0.0134 cbm
Weight of product 0.0008 kg
Weight of case 1.8 kg
EAN of inner box 4038917084321
EAN case 4038917003186

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

