



Product description

Order number	75.562.605
Product description	Container with screw cap VD, for virus diagnostics, (ØxH): 62 x 73 mm, transparent, graduated, screw cap, cap: light blue, cap assembled, non-sterile

Product characteristics

Graduation	yes
Type of cap	screw cap
Application	virus diagnostics

Size

Max. volume	100 ml
Diameter	62 mm
Diameter of opening	57 mm
Height of product	73 mm
Length excluding cap	72 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	light blue

Purity & certification

Product category	IVD product
CE certified	CE
Purity standard	non-sterile
Safety data sheet	nein
Batched	yes

Packaging

Minimum order qty.	200
Piece(s) / inner box	200
Piece(s) / case	200
Piece(s) / pallet	5600
Depth of case	588 mm
Width of case	398 mm
Height of case	255 mm
Case volume	0,0597 cbm
Weight of product	0,0153 kg
Weight of case	3,74 kg
EAN of inner box	4038917794053
EAN case	4038917794060

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

