



### **Product description**

Order number	75.562.605
Product description	Container with screw cap VD, for virus diagnostics, (ØxH): 62 x 73 mm, transparent, graduated, material: PP, screw cap, cap: light blue, cap assembled, non-sterile, 5 piece(s)/bag

### **Product characteristics**

Graduation	yes
Closure type	screw cap
Application	virus diagnostics

### **Size**

Max. volume	100 ml
Diameter	62 mm
Diameter of opening	57 mm
Height of product	73 mm
Length excluding cap	72 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	light blue

### **Purity & certification**

Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	non-sterile
Safety data sheet	no
Batched	yes

### **Packaging**

Minimum order qty.	200
Type of smallest subpackaging	bag
Piece(s) / inner box	200
Piece (s) / outer case	200
Piece(s) / pallet	5600
Depth of case	588 mm
Width of case	398 mm
Height of case	255 mm
Case volume	0.0597 cbm
Weight of product	0.0153 kg
Weight of case	3.74 kg
EAN of inner box	4038917794053
EAN case	4038917794060

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

