




**Product Specification**  
**Glass tubes, suitable for Coombs test,  
transparent**

**Page 1**



**Product description**

Order number 86.1509  
Product description Glass tubes, suitable for Coombs test, material: glass,  
(LxØ): 75 x 11.5 mm, transparent, 294 piece(s)/StackPack

**Product characteristics**

Application Coombs test  
Suitable for Coombs test  
Base shape round base

**Size**

Diameter 11,5 mm  
Length of product 75 mm

**Material & colours**

Product material glass  
Colour of product transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-10-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Glass tubes, suitable for Coombs test,**  
**transparent**

**Page 2**

**Purity & certification**

Product category	IVD product
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	2352
Type of inner packaging	StackPack
Piece(s) / inner box	2352
Piece(s) / case	2352
Piece(s) / pallet	56448
Case volume	0,0389 cbm
Weight of product	0,0038 kg
Weight of case	11,22 kg
EAN of inner box	4038917086783
EAN case	4038917072366

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-10-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**