



### **Product description**

Order number	95.1000.100
Product description	Tube segment opener, for bed side test cards, material: HD-PE, white, 200 piece(s)/case

### **Product characteristics**

Design	for bed side test cards
Application	Blood bank

### **Material & colours**

Product material	High Density Polyethylene (HD-PE)
Colour of product	white

### **Purity & certification**

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tube segment opener, for bed side test cards,**  
**HD-PE, white**

**Page 2**

**Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	case
Piece(s) / inner box	200
Piece (s) / outer case	2000
Piece(s) / pallet	242000
Depth of box	139 mm
Width of box	108 mm
Height of box	66 mm
Depth of case	338 mm
Width of case	230 mm
Height of case	160 mm
Case volume	0.0124 cbm
Weight of product	0.0018 kg
Weight of case	4.48 kg
EAN of inner box	4038917124607
EAN case	4038917028585

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-14

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**