



Product description

Order number	95.1000.100
Product description	Tube segment opener, for bed side test cards, material: HD-PE, white

Product characteristics

Design	for bed side test cards
Application	Blood bank

Material & colours

Product material	High Density Polyethylene (HD-PE)
Colour of product	white

Purity & certification

Product in accordance with MDD or IVDD	IVD product
CE certified	CE
Lot or serial number	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	2000
Type of inner packaging	case
Piece(s) / inner box	200
Piece(s) / case	2000
Piece(s) / pallet	242000
Depth of box	139 mm
Width of box	108 mm
Height of box	66 mm
Depth of case	338 mm
Width of case	228 mm
Height of case	160 mm
Case volume	0,0123 cbm
Weight of product	0,0018 kg
Weight of case	4,48 kg
EAN of inner box	4038917124607
EAN case	4038917028585

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

