




**Product Specification**  
**Tube segment opener, mounted on tube**  
**75 x 12 mm, HD-PE, white, tube: PS**

**Page 1**



**Product description**

Order number	95.1000.476
Product description	Tube segment opener, mounted on tube 75 x 12 mm, material: HD-PE, white, tube: PS, 50 piece(s)/StackPack

**Product characteristics**

Design	mounted on tube 75 x 12 mm
Application	Blood bank
Base shape	round base

**Size**

Diameter	12 mm
Length of product	75 mm

**Material & colours**

Product material	High Density Polyethylene (HD-PE)
Colour of product	white
Tube material	Polystyrene (PS)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-14

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**



## Product Specification

**Tube segment opener, mounted on tube  
75 x 12 mm, HD-PE, white, tube: PS**

**Page 2**

### Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

### Packaging

Minimum order qty.	400
Type of smallest subpackaging	StackPack
Piece(s) / inner box	50
Piece (s) / outer case	400
Piece(s) / pallet	57600
Depth of box	146 mm
Width of box	96 mm
Height of box	64 mm
Depth of case	323 mm
Width of case	213 mm
Height of case	145 mm
Case volume	0.01 cbm
Weight of product	0.0045 kg
Weight of case	2.22 kg
EAN of inner box	4038917124614
EAN case	4038917047975

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-14

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**