



Product description

Order number	95.1000.800
Product description	Tube segment opener, for glass tubes Ø 8,2 mm, material: HD-PE, white

Product characteristics

Application	Blood bank
Suitable for	glass tubes Ø 8,2 mm

Material & colours

Product material	High Density Polyethylene (HD-PE)
Colour of product	white

Purity & certification

Product in accordance with MDD or IVDD	IVD product
CE certified	CE
Lot or serial number	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	2000
Type of inner packaging	case
Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	216000
Depth of box	191 mm
Width of box	129 mm
Height of box	129 mm
Depth of case	398 mm
Width of case	268 mm
Height of case	150 mm
Case volume	0,016 cbm
Weight of product	0,0018 kg
Weight of case	4,38 kg
EAN of inner box	4038917124645
EAN case	4038917053273

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

