



Product Specification
Tube segment opener, for tubes with inner Ø 10 mm, HD-PE, white

Page 1



Product description

Order number 95.1000
Product description Tube segment opener, for tubes with inner Ø 10 mm, material: HD-PE, white, 500 piece(s)/case

Product characteristics

Application Blood bank
Suitable for tubes with inner Ø 10 mm

Size

Diameter 10 mm

Material & colours

Product material High Density Polyethylene (HD-PE)
Colour of product white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-04

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Tube segment opener, for tubes with inner Ø
10 mm, HD-PE, white

Page 2

Purity & certification

Product category	In-vitro diagnostic, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	case
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	216000
Depth of box	191 mm
Width of box	129 mm
Height of box	129 mm
Depth of case	398 mm
Width of case	268 mm
Height of case	150 mm
Case volume	0.016 cbm
Weight of product	0.0015 kg
Weight of case	3.8 kg
EAN of inner box	4038917084949
EAN case	4038917053266

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-04

This document was prepared by EDP support and is valid without signature.



SARSTEDT