




**Product Specification**  
**Absorbent liner, suitable for Mailing container 114  
x 44 mm, (LxW): 135 x 95 mm**

**Page 1**



**Product description**

Order number	78.669
Product description	Absorbent liner, suitable for Mailing container 114 x 44 mm, (LxW): 135 x 95 mm, 100 piece(s)/case

**Product characteristics**

Application	Mailing of samples, Multi-Safe disposel boxes
Suitable for	Mailing container 114 x 44 mm

**Size**

Width of product	95 mm
Length of product	135 mm

**Material & colours**

Product material	cellulose
------------------	-----------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Absorbent liner, suitable for Mailing container 114**  
**x 44 mm, (LxW): 135 x 95 mm**

**Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

**Packaging**

Minimum order qty.                      1000  
Type of smallest                              case  
subpackaging  
Piece(s) / inner box                      100  
Piece (s) / outer case                      1000  
Piece(s) / pallet                              56000  
Depth of case                              385 mm  
Width of case                                385 mm  
Height of case                               160 mm  
Case volume                                0.0237 cbm  
Weight of product                        0.0026 kg  
Weight of case                              3.04 kg  
EAN of inner box                            4038917083249  
EAN case                                    4038917002196

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

