




**Product Specification**  
**Absorbent liner, suitable for cold transport container, (LxW): 80 x 55 mm**

**Page 1**



**Product description**

|                     |  |
|---------------------|--|
| Order number        | 78.671   |
| Product description | Absorbent liner, suitable for cold transport container, (LxW): 80 x 55 mm, 100 piece(s)/case |

**Product characteristics**

|              |   |
|--------------|---|
| Application  | Mailing of samples, Multi-Safe disposel boxes |
| Suitable for | cold transport container                      |

**Size**

|                   |       |
|-------------------|-------|
| Width of product  | 55 mm |
| Length of product | 80 mm |

**Material & colours**

|                  |           |
|------------------|-----------|
| Product material | cellulose |
|------------------|-----------|

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Absorbent liner, suitable for cold transport container, (LxW): 80 x 55 mm**

**Page 2**

**Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | no                      |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 2500          |
| Type of smallest subpackaging | case          |
| Piece(s) / inner box          | 100           |
| Piece (s) / outer case        | 2500          |
| Piece(s) / pallet             | 270000        |
| Depth of case                 | 383 mm        |
| Width of case                 | 193 mm        |
| Height of case                | 190 mm        |
| Case volume                   | 0.014 cbm     |
| Weight of product             | 0.0009 kg     |
| Weight of case                | 2.57 kg       |
| EAN of inner box              | 4038917083263 |
| EAN case                      | 4038917079860 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-08

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**