



### **Product description**

Order number 93.1650.010  
Product description Insertion plate, Version 1, material: PP, 1 piece(s)/case

### **Product characteristics**

Suitable for Safety Tray (93.1650) & e.g. 1x Multi-Safe twin plus 500 (77.3884.005/.105), 1x rack with handle (93.851 - 93.851.1xx), 2x dispenser boxes (95.813.001), 3x multi-purpose cup (75.1355), 1x disinfection spray

### **Size**

Width of product 272 mm  
Length of product 420 mm

### **Material & colours**

Product material Polypropylene (PP)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1440
Case volume	0.0009 cbm
Weight of product	0.152 kg
Weight of case	0.31 kg
EAN of inner box	4038917132695
EAN case	4038917076715

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

