



### **Product description**

Order number 95.813.001  
Product description Dispenser box, material: PS, lid assembled

### **Product characteristics**

Base shape flat base

### **Size**

Diameter of base 81 mm  
Diameter of lid 83 mm  
Diameter of opening 81 mm  
Length of product 110 mm

### **Material & colours**

Product material Polystyrene (PS)  
Colour of product white  
Closure material Polystyrene (PS)  
Colour of cap white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category                   no medical device | IVD  
Batched                               no

**Packaging**

Minimum order qty.               35  
Piece(s) / inner box               35  
Piece (s) / outer case             35  
Piece(s) / pallet                  1260  
Depth of case                    398 mm  
Width of case                     398 mm  
Height of case                    300 mm  
Case volume                      0.0475 cbm  
Weight of product                0.0623 kg  
Weight of case                    2.77 kg  
EAN of inner box                 4038917121644  
EAN case                         4038917074797

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

