



### **Product description**

Order number 95.813.001  
Product description Dispenser box material: PS

### **Product characteristics**

Base shape flat base

### **Size**

Diameter of base 81 mm  
Diameter of lid 83 mm  
Diameter of opening 81 mm  
Length of product 110 mm

### **Material & colours**

Product material Polystyrene (PS)  
Colour of product white  
Cap material Polystyrene (PS)  
Colour of cap white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Product category                   no medical device | IVD  
Batched                               no

### **Packaging**

Minimum order qty.               35  
Piece(s) / inner box               35  
Piece(s) / case                   35  
Piece(s) / pallet                 1260  
Depth of case                   398 mm  
Width of case                   398 mm  
Height of case                  300 mm  
Case volume                   0,0475 cbm  
Weight of product               0,0623 kg  
Weight of case                 2,77 kg  
EAN of inner box               4038917121644  
EAN case                       4038917074797

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

