



Product description

Order number	95.963.001
Product description	Wall holder, for Multi-Safe euroMatic® , material: stainless steel

Product characteristics

Suitable for	Multi-Safe euroMatic®
--------------	-----------------------

Size

Width of product	161 mm
Length of product	222 mm

Material & colours

Product material	stainless steel
Colour of product	silver/white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest case
subpackaging
Piece(s) / inner box 1
Piece (s) / outer case 20
Piece(s) / pallet 2880
Depth of case 323 mm
Width of case 213 mm
Height of case 145 mm
Case volume 0.01 cbm
Weight of product 0.246 kg
Weight of case 5.24 kg
EAN of inner box 4038917233620
EAN case 4038917233613

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

