




Product Specification

Page 1

Fastening plate for Dräger rail system and combinable with euroMatic wall holder® (95.963.001), vario 2000 wall holder (95.963.020) and twin plus wall holder (95.963.007)



Product description

Order number	95.963.005
Product description	Fastening plate, for Dräger rail system and combinable with euroMatic wall holder® (95.963.001), vario 2000 wall holder (95.963.020) and twin plus wall holder (95.963.007), material: stainless steel, (LxWxH): 130 x 15 x 79 mm

Product characteristics

Application	For attachment to system bed rails in combination with wall holder euroMatic® (95.963.001), wall holder vario 2000 (95.963.020) and wall holder twin plus (95.963.007)
Suitable for	Dräger rail system and combinable with euroMatic wall holder® (95.963.001), vario 2000 wall holder (95.963.020) and twin plus wall holder (95.963.007)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

Page 2

Fastening plate for Dräger rail system and combinable with euroMatic wall holder® (95.963.001), vario 2000 wall holder (95.963.020) and twin plus wall holder (95.963.007)

Size

Width of product	15 mm
Height of product	79 mm
Length of product	130 mm

Material & colours

Product material	stainless steel
Colour of product	grey

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	20
Piece(s) / pallet	6480
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.1752 kg
Weight of case	3.68 kg
EAN of inner box	4038917364096
EAN case	4038917364089

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT