




## Product Specification

**Fastening plate for Dräger rail system and combinable with wall holder euroMatic® (95.963.001), wall holder vario 2000 (95.963.020) and wall holder twin plus (95.963.007)**

**Page 1**



### Product description

Order number 95.963.005  
Product description Fastening plate, for Dräger rail system and combinable with wall holder euroMatic® (95.963.001), wall holder vario 2000 (95.963.020) and wall holder twin plus (95.963.007), material: stainless steel, (LxWxH): 130 x 15 x 79 mm

### Product characteristics

Application For attachment to system bed rails in combination with wall holder euroMatic® (95.963.001), wall holder vario 2000 (95.963.020) and wall holder twin plus (95.963.007)  
Suitable for Dräger rail system and combinable with wall holder euroMatic® (95.963.001), wall holder vario 2000 (95.963.020) and wall holder twin plus (95.963.007)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2023-11-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Fastening plate for Dräger rail system and combinable with wall holder euroMatic® (95.963.001), wall holder vario 2000 (95.963.020) and wall holder twin plus (95.963.007)**

Page 2

### Size

Width of product	15 mm
Height of product	79 mm
Length of product	130 mm

### Material & colours

Product material	stainless steel
Colour of product	grey

### Purity & certification

Product category	no medical device   IVD
Batched	no

### Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	20
Piece(s) / pallet	6480
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.1752 kg
Weight of case	3.68 kg
EAN of inner box	4038917364096
EAN case	4038917364089

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2023-11-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**