



## **Product description**

Order number	77.1515
Product description	Disposal container, 2 l, transparent, label/print: yellow/black, material: PP, opening: round, height: 245 mm, flat cap, Ø of cap: 82 mm, cap: yellow, 30 piece(s)/carton

## **Product characteristics**

Colour of print/label	yellow/black
Cap	flat cap
Closure type	screw cap
Shape of opening	round

## **Size**

Max. volume	2 l
Diameter	116 mm
Diameter of base	105 mm
Diameter of lid	82 mm
Diameter of opening	71 mm
Height of product	245 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Colour of cap	yellow

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	30
Type of smallest subpackaging	case
Piece(s) / inner box	30
Piece (s) / outer case	30
Piece(s) / pallet	360
Depth of case	608 mm
Width of case	368 mm
Height of case	515 mm
Case volume	0.1152 cbm
Weight of product	0.125 kg
Weight of case	4.5 kg
EAN of inner box	4038917086813
EAN case	4038917005319

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

