



Product description

Order number	77.3884.105
Product description	Multi-Safe twin plus 500, 500 ml, Biohazard, material: PP, application e.g.: disposal of pointed and sharp objects, opening: oval, height: 162 mm, 126 piece(s)/carton

Product characteristics

Label/ Print	with barcode label
Colour of print/label	yellow/red
Example of use	disposal of pointed and sharp objects
Shape of opening	oval
Removal aid - turning aid	Removal and rotation aid

Size

Max. volume	500 ml
Diameter	89 mm
Diameter of base	57 mm
Diameter of lid	89 mm
Height of product	162 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polypropylene (PP)
Colour of product	yellow
Colour of cap	red
Colour of lid	red

Purity & certification

Product in accordance with MDD or IVDD	no medical device IVD
Lot or serial number	no

Packaging

Minimum order qty.	126
Piece(s) / inner box	126
Piece(s) / case	126
Piece(s) / pallet	3024
Case volume	0,0706 cbm
Weight of product	0,0753 kg
Weight of case	10,24 kg
EAN of inner box	4038917542463
EAN case	4038917542456

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

