



Product description

Order number	77.3884.130
Product description	Disposal container, Multi-Safe twin plus, 3,000 ml, yellow, with plastic label, label/print: yellow/red, biohazard labelling, material: PP, application e.g.: disposal of pointed and sharp objects, opening: oval, height: 233 mm, lid (LxW): 179 x 176 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 62 piece(s)/carton

Product characteristics

Label/ Print	with plastic label
Colour of print/label	yellow/red
Example of use	disposal of pointed and sharp objects
Shape of opening	oval
Removal aid - turning aid	Removal and rotation aid

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.



Size

Max. volume	3000 ml
Diameter of base	112 mm
Height of product	233 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	yellow
Colour of cap	red
Colour of lid	red

Purity & certification

Satisfies the requirement	TRBA 250, UN 3291, DIN EN ISO 23907
Product category	not regulated
CE certified	kein CE-Kennzeichen
Batched	no

Packaging

Minimum order qty.	62
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	62
Piece(s) / pallet	744
Depth of box	180 mm
Width of box	180 mm
Height of box	238 mm
Depth of case	800 mm
Width of case	400 mm
Height of case	380 mm
Case volume	0.1216 cbm
Weight of product	0.241 kg
Weight of case	16.53 kg
EAN of inner box	4038917903523
EAN case	4038917903530

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

