



### **Product description**

Order number	77.3890.002
Product description	Disposal container, Multi-Safe sani, 200 ml, yellow, material: PP, opening: oval, height: 131 mm, flat cap, lid (LxW): 75 x 47 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 110 piece(s)/carton
UN identification	yes

### **Product characteristics**

Cap	flat cap
Shape of opening	oval
Removal aid - turning aid	Removal aid

### **Size**

Max. volume	200 ml
Diameter of base	57 mm x 32 mm
Diameter of opening	33 mm x 23 mm
Height of product	131 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polypropylene (PP)
Colour of product	yellow
Colour of cap	red
Colour of lid	red

**Purity & certification**

Satisfies the requirement	TRBA 250, UN 3291, DIN EN ISO 23907
Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	110
Type of smallest subpackaging	case
Piece(s) / inner box	110
Piece (s) / outer case	110
Piece(s) / pallet	3300
Depth of case	397 mm
Width of case	325 mm
Height of case	397 mm
Case volume	0.0512 cbm
Weight of product	0.0483 kg
Weight of case	6 kg
EAN of inner box	4038917330916
EAN case	4038917330909

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

