



Product description

Order number	77.3890.202
Product description	Disposal container, Multi-Safe sani REC, 200 ml, grey, material: recycled PP, opening: oval, height: 131 mm, lid (LxW): 75 x 47 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 110 piece(s)/carton

Product characteristics

Shape of opening	oval
Removal aid - turning aid	Removal aid

Size

Max. volume	200 ml
Diameter of base	57 mm x 32 mm
Diameter of opening	33 mm x 23 mm
Height of product	131 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	recycled Polypropylene (recycled PP)
Colour of product	grey
Colour of cap	red
Colour of lid	red

Purity & certification

Satisfies the requirement	TRBA 250, UN 3291, DIN EN ISO 23907
Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Minimum order qty.	110
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	110
Piece(s) / pallet	3300
Depth of case	397 mm
Width of case	397 mm
Height of case	325 mm
Case volume	0.0512 cbm
Weight of product	0.0483 kg
Weight of case	6 kg
EAN of inner box	4066916004484
EAN case	4066916004491

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

