



Product description

| | |
|---------------------|--|
| Order number | 77.3893.120 |
| Product description | Disposal container, Multi-Safe vario, 2,000 ml, yellow, square, with plastic label, label/print: yellow/red, biohazard labeling, material: PP, application e.g.: disposal of pointed and sharp objects, opening: oval, height: 140 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 50 piece(s)/carton |

Product characteristics

| | |
|---------------------------|---------------------------------------|
| Label/ Print | with plastic label |
| Colour of print/label | yellow/red |
| Width of label | 160 mm |
| Height of label | 70 mm |
| Example of use | disposal of pointed and sharp objects |
| Absorbent liner | no |
| Shape of opening | oval |
| Removal aid - turning aid | Removal and rotation aid |
| Liquid-tight | no |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|---------------------|---------------|
| Max. volume | 2000 ml |
| Diameter of opening | 95 mm x 60 mm |
| Width of product | 95 mm |
| Height of product | 140 mm |
| Length of product | 240 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | yellow |
| Closure material | Polypropylene (PP) |
| Colour of cap | red |

Purity & certification

| | |
|---------------------------|-------------------------------------|
| Satisfies the requirement | TRBA 250, UN 3291, DIN EN ISO 23907 |
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 50 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 50 |
| Piece(s) / pallet | 1000 |
| Depth of case | 400 mm |
| Width of case | 365 mm |
| Height of case | 600 mm |
| Case volume | 0.0876 cbm |
| Weight of product | 0.2272 kg |
| Weight of case | 12.4 kg |
| EAN of inner box | 4038917532419 |
| EAN case | 4038917540230 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

