



Product description

Order number	77.3898.011
Product description	Disposal container, Multi-Safe medi, 11 l, yellow, round, material: PP, opening: round, height: 220 mm, flat cap, Ø of cap: 290 mm, cap: red, TRBA 250, DIN EN ISO 23907, 250 piece(s)/carton
UN identification	no

Product characteristics

Cap	flat cap
Shape of opening	round
Removal aid - turning aid	Removal aid

Size

Max. volume	11 l
Diameter of base	257 mm
Diameter of lid	290 mm
Diameter of opening	88 mm
Height of product	220 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	yellow
Colour of cap	red
Colour of lid	yellow

Purity & certification

Satisfies the requirement	TRBA 250, DIN EN ISO 23907
Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	250
Type of smallest subpackaging	pallet
Piece(s) / inner box	250
Piece (s) / outer case	250
Piece(s) / pallet	250
Depth of case	792 mm
Width of case	1090 mm
Height of case	1192 mm
Case volume	1.029 cbm
Weight of product	0.5121 kg
Weight of case	137.97 kg
EAN of inner box	4038917409483
EAN case	4038917409476

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

