



## **Product description**

Order number	77.3898.206
Product description	Disposal container, Multi-Safe medi, 6 l, yellow, round, with plastic label, label/print: black, material: PP, with assembled bracket, opening: round, height: 195 mm, flat cap, Ø of cap: 230 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 21 piece(s)/carton
UN identification	yes

## **Product characteristics**

Label/ Print	with plastic label
Colour of print/label	black
Cap	flat cap
Design	with assembled bracket
Shape of opening	round
Removal aid - turning aid	Removal and rotation aid

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

Max. volume	6 l
Diameter of base	190 mm
Diameter of lid	230 mm
Diameter of opening	88 mm
Height of product	195 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	yellow
Colour of cap	red
Colour of lid	yellow

### **Purity & certification**

Satisfies the requirement	TRBA 250, UN 3291, DIN EN ISO 23907
Product category	not regulated
CE certified	no CE mark
Batched	no

### **Packaging**

Minimum order qty.	21
Type of smallest subpackaging	case
Piece(s) / inner box	21
Piece (s) / outer case	21
Piece(s) / pallet	336
Depth of box	395 mm
Width of box	430 mm
Height of box	585 mm
Depth of case	585 mm
Width of case	395 mm
Height of case	430 mm
Case volume	0.0994 cbm
Weight of product	0.3472 kg
Weight of case	8 kg
EAN of inner box	4038917543453
EAN case	4038917543460

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

