



Product description

Order number	77.3898.511
Product description	Disposal container, Multi-Safe medi, 11 l, yellow, round, material: PP, opening: round, height: 220 mm, flat cap, Ø of cap: 290 mm, cap: red, TRBA 250, DIN EN ISO 23907, 10 piece(s)/carton
UN identification	no

Product characteristics

Cap	flat cap
Shape of opening	round
Removal aid - turning aid	Removal aid

Size

Max. volume	11 l
Diameter of base	257 mm
Diameter of lid	290 mm
Diameter of opening	88 mm
Height of product	220 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	yellow
Colour of cap	red
Colour of lid	yellow

Purity & certification

Satisfies the requirement	TRBA 250, DIN EN ISO 23907
Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Minimum order qty.	10
Type of smallest subpackaging	case
Piece(s) / inner box	10
Piece (s) / outer case	10
Piece(s) / pallet	160
Depth of box	395 mm
Width of box	430 mm
Height of box	585 mm
Depth of case	585 mm
Width of case	395 mm
Height of case	430 mm
Case volume	0.0994 cbm
Weight of product	0.637 kg
Weight of case	6.87 kg
EAN of inner box	4038917409544
EAN case	4038917409537

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

