



### **Product description**

Order number	77.3899.035
Product description	Disposal container, Multi-Safe steri, 35 l, yellow, square, material: PP, autoclavable, opening: square, height: 330 mm, flat cap, Cover with seal, lid (LxW): 393 x 393 mm, cap: red, UN 3291, 150 piece(s)/carton
UN identification	yes

### **Product characteristics**

Cap	flat cap
Shape of opening	square
Removal aid - turning aid	not available
Autoclavable (max. °C)	134 °C

### **Size**

Max. volume	35 l
Diameter of base	340 mm x 340 mm
Diameter of opening	370 mm x 370 mm
Height of product	330 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	yellow
Colour of cap	red
Colour of lid	yellow

### **Purity & certification**

Satisfies the requirement	UN 3291
Product category	not regulated
CE certified	no CE mark
Batched	no

### **Packaging**

Minimum order qty.	150
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	150
Piece(s) / pallet	150
Depth of case	1200 mm
Width of case	800 mm
Height of case	2430 mm
Case volume	2.304 cbm
Weight of product	1.1945 kg
Weight of case	179.25 kg
EAN of inner box	4038917333498
EAN case	4038917333481

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

