



Product description

Order number	77.3899.060
Product description	Disposal container, Multi-Safe steri, 60 l, yellow, square, material: PP, autoclavable, opening: square, height: 520 mm, flat cap, lid (LxW): 393 x 393 mm, cap: red, UN 3291, 120 piece(s)/carton
UN identification	yes

Product characteristics

Cap	flat cap
Shape of opening	square
Removal aid - turning aid	not available
Autoclavable (max. °C)	134 °C

Size

Max. volume	60 l
Diameter of base	340 mm x 340 mm
Diameter of opening	370 mm x 370 mm
Height of product	520 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	yellow
Colour of cap	red
Colour of lid	yellow

Purity & certification

Satisfies the requirement	UN 3291
Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Minimum order qty.	120
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	120
Piece(s) / pallet	120
Depth of case	1200 mm
Width of case	800 mm
Height of case	2260 mm
Case volume	2.304 cbm
Weight of product	1.8345 kg
Weight of case	220.2 kg
EAN of inner box	4038917333559
EAN case	4038917333542

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

