



Product description

Order number	51.1534.500
Product description	Salivette® Cortisol, with synthetic swab, cap: blue, with paper label, label/print: white/blue, for cortisol determination from saliva, 100 piece(s)/bag

Product characteristics

Label/ Print	with paper label
Colour of print/label	white/blue
Closure type	push cap
Design	with synthetic swab
Application	cortisol determination from saliva
Recommended centrifugation	1000 x g
Centrifugation time	2 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Diameter	16.8 mm
Length including cap	97 mm
Length excluding cap	84 mm

Material & colours

Product material	Polypropylene (PP)
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	blue

Purity & certification

Product category	In-vitro diagnostic, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	500
Type of inner packaging	bag
Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	28000
Depth of box	226 mm
Width of box	101 mm
Height of box	119 mm
Depth of case	533 mm
Width of case	263 mm
Height of case	140 mm
Case volume	0.0196 cbm
Weight of product	0.0069 kg
Weight of case	4.09 kg
EAN of inner box	4038917159364
EAN case	4038917159357

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

