



## Product description

|                     |  |
|---------------------|--|
| Order number        | 51.1534.500  |
| Product description | Salivette® Cortisol, with synthetic swab, cap: blue, with paper label, label/print: white/blue, for cortisol determination from saliva, 100 piece(s)/bag |

## Product characteristics

|                            |                                    |
|----------------------------|------------------------------------|
| Label/ Print               | with paper label                   |
| Colour of print/label      | white/blue                         |
| Closure type               | push cap                           |
| Design                     | with synthetic swab                |
| Application                | cortisol determination from saliva |
| Recommended centrifugation | 1000 x g                           |
| Centrifugation time        | 2 Min                              |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

|                      |         |
|----------------------|---------|
| Diameter             | 16.8 mm |
| Length including cap | 97 mm   |
| Length excluding cap | 84 mm   |

### **Material & colours**

|                  |                                  |
|------------------|----------------------------------|
| Product material | Polypropylene (PP)               |
| Closure material | Low Density Polyethylene (LD-PE) |
| Colour of cap    | blue                             |

### **Purity & certification**

|                  |                               |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified     | CE                            |
| Batched          | yes                           |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 500           |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 100           |
| Piece (s) / outer case        | 500           |
| Piece(s) / pallet             | 28000         |
| Depth of box                  | 226 mm        |
| Width of box                  | 101 mm        |
| Height of box                 | 119 mm        |
| Depth of case                 | 533 mm        |
| Width of case                 | 263 mm        |
| Height of case                | 140 mm        |
| Case volume                   | 0.0196 cbm    |
| Weight of product             | 0.0069 kg     |
| Weight of case                | 4.09 kg       |
| EAN of inner box              | 4038917159364 |
| EAN case                      | 4038917159357 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

