



Product Specification
Salivette® , with cotton swab, cap: white, with
paper label

Page 1



Product description

Order number	51.1534
Product description	Salivette® , with cotton swab, cap: white, with paper label, label/print: white/black, for Saliva examination, 100 piece(s)/bag

Product characteristics

Label/ Print	with paper label
Colour of print/label	white/black
Closure type	push cap
Design	with cotton swab
Application	Saliva examination
Recommended centrifugation	1000 x g
Centrifugation time	2 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Salivette® , with cotton swab, cap: white, with paper label

Page 2

Size

Diameter	16.8 mm
Length including cap	97 mm
Length excluding cap	84 mm

Material & colours

Product material	Polypropylene (PP)
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	white

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	28000
Depth of case	598 mm
Width of case	398 mm
Height of case	125 mm
Case volume	0.0298 cbm
Weight of product	0.007 kg
Weight of case	4.09 kg
EAN of inner box	4038917086905
EAN case	4038917005432

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

