



Product description

Order number	71.9921.316
Product description	Bucket, up to 1 l, (ØxH): 122 x 145 mm, white, cap enclosed, for collection and storage of diagnostic specimens, 24 piece(s)/case

Product characteristics

Application	collection and storage of diagnostic specimens
Autoclavable (max. °C)	121 °C

Size

Max. volume	1 l
Diameter	122 mm
Height of product	145 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Latex-free	yes
Pyrogen-free	no
Safety data sheet	nein

Packaging

Minimum order qty.	24
Type of inner packaging	case
Piece(s) / inner box	24
Piece(s) / case	24
Piece(s) / pallet	1152
Depth of case	428 mm
Width of case	358 mm
Height of case	234 mm
Case volume	0,0357 cbm
Weight of case	2,7 kg
EAN of inner box	4038917150620
EAN case	4038917038195

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

