



### **Product description**

Order number	71.9921.317
Product description	Bucket, up to 2 l, (ØxH): 154 x 180 mm, white, cap enclosed, for collection and storage of diagnostic specimens, 24 piece(s)/case

### **Product characteristics**

Application	collection and storage of diagnostic specimens
Autoclavable (max. °C)	121 °C

### **Size**

Max. volume	2 l
Diameter	154 mm
Height of product	180 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Latex-free	yes
Pyrogen-free	no
Safety data sheet	no

### **Packaging**

Minimum order qty.	24
Type of smallest subpackaging	case
Piece(s) / inner box	24
Piece (s) / outer case	24
Piece(s) / pallet	960
Depth of case	604 mm
Width of case	264 mm
Height of case	373 mm
Case volume	0.06 cbm
Weight of case	4.6 kg
EAN of inner box	4038917150637
EAN case	4038917038201

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

