



### **Product description**

Order number	71.9921.318
Product description	Bucket, up to 5 l, (ØxH): 189 x 200 mm, white, cap enclosed, for collection and storage of diagnostic specimens, 12 piece(s)/case

### **Product characteristics**

Application	collection and storage of diagnostic specimens
Autoclavable (max. °C)	121 °C

### **Size**

Max. volume	5 l
Diameter	189 mm
Height of product	200 mm
Length of product	223,5 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	white

### **Purity & certification**

Latex-free	yes
Pyrogen-free	no
Safety data sheet	nein

### **Packaging**

Minimum order qty.	12
Type of inner packaging	case
Piece(s) / inner box	12
Piece(s) / case	12
Piece(s) / pallet	360
Depth of case	548 mm
Width of case	373 mm
Height of case	408 mm
Case volume	0,0834 cbm
Weight of case	4,2 kg
EAN of inner box	4038917150644
EAN case	4038917398381

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

