



Product description

Order number	71.9921.318
Product description	Bucket, up to 5 l, (ØxH): 189 x 200 mm, white, cap enclosed, for collection and storage of diagnostic specimens, 12 piece(s)/case

Product characteristics

Application	collection and storage of diagnostic specimens
Autoclavable (max. °C)	121 °C

Size

Max. volume	5 l
Diameter	189 mm
Height of product	200 mm
Length of product	223.5 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product white

Purity & certification

Latex-free yes
Pyrogen-free no
Safety data sheet no

Packaging

Minimum order qty. 12
Type of smallest subpackaging case
Piece(s) / inner box 12
Piece (s) / outer case 12
Piece(s) / pallet 360
Depth of case 548 mm
Width of case 373 mm
Height of case 408 mm
Case volume 0.0834 cbm
Weight of case 4.2 kg
EAN of inner box 4038917150644
EAN case 4038917398381

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

